

PAL-650 Cleaning & Sterilization

CLEAN (CAP ON)

- ▶ Connect CAP-600E or Cable to keep liquid out.
- ▶ Scrub with warm water, a soft brush and pH neutral detergent, then rinse.*
- ▶ *Do not immerse. Do not put nose of handpiece into water stream.*



STERILIZE (CAP OFF)

- ▶ Remove CAP-600E or Cable to allow steam to penetrate handpiece.
- ▶ Full dry cycle is required to prevent corrosion and to ensure proper performance.
- ▶ Cool slowly at room temperature. **Never run while warm.**
- ▶ Inspect carefully to ensure all moisture is completely removed.**



*Exposure Time
(Full Cycle)*



Temperature



Dry Time



Inspect

Dynamic-Air-Removal

4 minutes

270° F
(132° C)

20 minutes



Dry Instrument - Okay to Use