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ENGLISH

ULTRATINE™ Forehead Device Instructions for Use

Available Models		
Part Number	Time Length	Part Name
23202	3.0 mm	ULTRATINE Forehead 3.0
23102	3.5 mm	ULTRATINE Forehead 3.5

DESCRIPTION
The ULTRATINE™ Forehead Device is a bioabsorbable fixation implant. This device is SINGLE USE ONLY and is supplied sterile. An Instrument Kit specific for this device is packaged and supplied separately.

INDICATIONS
ULTRATINE Forehead is intended for use in subperiosteal browplasty surgery to fixate the sub-dermis to the cranial bone and to maintain fixation until healing has occurred.

CONTRAINDIKATIONER
1. **Patients appearing to have very thin bones which might imply an inadequate cranial thickness for the ULTRATINE Forehead Device bone post which extends to a depth of 3.75 mm. The drilled hole extends to a depth of 3.95 mm +0.05, -0.0 (max 4.00 mm).**
2. Situations where internal fixation is otherwise contraindicated (e.g. infection).
3. Thin, atrophic scalp.
4. Any known allergy or foreign-body sensitivities to plastic biomaterial, such as PLA/PGA.

LABELLING, PACKAGING AND STERILIZATION
ULTRATINE Forehead Device is sterilized by gamma irradiation. **DO NOT** sterilize this device. Implant should be accepted only if the factory packaging and labeling arrive intact. Do not use if packaging shows evidence of punctures, tampering, water contamination or other damage.

STORAGE INSTRUCTIONS
1. Store at room temperature (15° to 24° C or 60° to 75° F) in a dry place out of direct sunlight. Do not use beyond the expiration date listed on the label. Do not use if the temperature indicator dot has turned black, the indicator changes color if the device is subjected to temperatures that could damage it, even if the exposure is of short duration, as might occur during transportation.

INSTRUCTIONS FOR USE (for endoscopic or open techniques)
1. Prior to patient usage, verify separately that the scalp thickness is appropriate for the ULTRATINE Forehead implants.
2. Make incisions large enough to allow insertion of the ULTRATINE Forehead device without bending during placement, approximately 20 mm.
3. Perform the dissection in the usual fashion, making sure to achieve adequate release.
4. To minimize the risk of creating a full thickness cranial hole with the Forehead drill bit, verify that the hole position is medial to the temporal fusion line and anterior to the coronal suture (**Figure 1**). Avoid the sagittal sinus at midline.
5. Plan the cranial hole placement such that the implant will ultimately lie under hair bearing tissue. Typically, for sagittal, paramedian incisions, the implant lies somewhat lateral and/or anterior to the incision. For coronal or transverse incisions, the implant lies anterior to the incision.
6. Prepare the cranial bone using an ENDOTINE Forehead™ drill bit (**Figure 2**) applied perpendicular to the bone. **After drilling, check that no residual loose bone material is present in the drill hole.**
7. Remove the ULTRATINE Forehead Device from the package foam using the ENDOTINE Forehead™ Insertion Tool. Properly align the insertion tool with the ULTRATINE Forehead Device to avoid damage (**Figure 3**).
8. Utilizing the insertion tool, align the post of the device perpendicular to and within the center of the hole before pushing straight down (**Figure 4**).
9. With moderate pressure, insert the ULTRATINE Forehead Device into the prepared site until the device platform is resting on the cranium (**Figure 5**).

10. Release the ULTRATINE Forehead Device from the Insertion Tool. **The device cannot be removed from the hole and reused.**
11. Perform the brow elevation and apply manual pressure to the scalp to ensure penetration of the tissue by the device (**Figure 6**).
12. If a drain is used, manipulate the subcutaneous portion such that it is maximal distance from the implant.
13. Use caution during skin closure to avoid elevating the scalp tissue off the implant.
14. A gentle pressure dressing may be appropriate upon completion of the case.

PRECAUTIONS
1. The implant ideally lies under hair bearing scalp at the conclusion of the procedure. Placement under forehead or alopectic skin may lead to unacceptable visibility of the device prior to absorption.
2. The patient is to be informed of possible adverse effects including palpability of the device prior to expiration, discomfort, and surgical risks. Tenderness and visible fullness over the device for several months has been noted, especially in patients with thin scalps. These symptoms and/or wound problems may lead to a decision to surgically remove the device prior to full resorption.
3. Local cyst or sinus formation may occur and may be treated by aspiration or debridement.
4. Discard open, unused devices.
5. This device is for SINGLE USE ONLY and **MUST NEVER BE RE-USED.** Restertization may cause loss of function or failure of the device. Risks of re-use may include patient infection due to lack of sterility.
6. Visually inspect the insertion tool and drill bits for erosion, cracks or damage after each sterilization and prior to use. Do not attempt to repair or sharpen instruments.
7. Incomplete insertion of the device in the cranial bone may result in inadequate anchoring.
8. Patients should be instructed to avoid contact sports or trauma to the head during healing.
9. Contour irregularities of the skin may be noticeable overlying the ULTRATINE Forehead Device, especially in non hair-bearing areas.

ADVERSE EFFECTS
Material sensitivity/allergic reactions in patients following surgery should be reported. Implantation of foreign materials in tissues can result in histological reactions. In rare cases, the device may result in granuloma or become encapsulated. If granuloma or encapsulation is suspected, prompt treatment is recommended.
Caution: Federal Law (USA) restricts this device to sale, distribution, or use by or on the order of a physician.

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