

MicroAire Surgical Instruments, LLC
3590 Grand Forks Boulevard
Charlottesville, Virginia 22911 USA

Phone: 434 975 8000
Fax: 434 975 4131
www.microaire.com



MRI Compatibility

October 29, 2015

To Whom It May Concern,

MicroAire Surgical Instruments does not recommend a patient undergo an MRI procedure while a Kirschner Wire or Steinmann pin is fixated in their body, as these wires and pins have a potential to move during an MRI procedure. Please consult with your healthcare provider prior to undergoing an MRI procedure.

If you have any questions or concerns, please contact MicroAire Quality Assurance & Regulatory Affairs Department at the phone or fax number listed above.

Kind Regards,

A handwritten signature in black ink, appearing to read 'Glenn Gerstenfeld', is positioned above the printed name.

Glenn Gerstenfeld
Vice President of Quality Assurance and Regulatory Affairs