

MicroAire Surgical Instruments, LLC
3590 Grand Forks Boulevard
Charlottesville, Virginia 22911 USA

Phone: 434 975 8000
Fax: 434 975 4131
www.microaire.com



Nickel Content

October 29, 2015

To Whom It May Concern,

MicroAire Surgical Instruments is aware of concerns regarding the nickel content of Kirschner Wires and Steinmann Pins. We recommend that a patient with a potential metal sensitivity be seen by a dermatologist or allergist and undergo appropriate testing before having surgery.

Nickel Composition (weight percent): 13.00 to 15.00

If you have any questions or concerns, please contact MicroAire Quality Assurance & Regulatory Affairs Department at the phone or fax number listed above.

Kind Regards,

A handwritten signature in black ink, appearing to read 'Glenn Gerstenfeld', is written above the printed name.

Glenn Gerstenfeld

Vice President of Quality Assurance and Regulatory Affairs