

POWERFUL, SMOOTH, FAST, AND QUIET



MICROAIRE 5130 High-Speed Drill

- ✓ Connects to both electric or pneumatic Series-5000 Small-Bone Handpieces.
- ✓ MicroAire® SmartLoop technology maintains constant 83K RPM, even under load.
- ✓ Works with hand-lever throttle or 5401 footswitch.
- ✓ Comes with an extra medium bur guard, cleaning brush, and replacement O-rings.
- ✓ Long bur guard also available.
- ✓ Uses MicroAire® ZB-series burs (excluding extra-long burs).

ACCESSORIES

5130-BG-M	Bur Guard, Medium
5130-BG-L	Bur Guard, Long
5130-BR	Cleaning brush
5130-SP	Spare Parts Kit (O-rings)

Contact
to schedule a product evaluation.

microaire.com

© MicroAire Surgical Instruments LLC. All rights reserved.

MicroAire Surgical Instruments
3590 Grand Forks Boulevard, Charlottesville, VA 22911 USA

LIT-5130 Rev A (June 2021)

MICROAIRE®
For Surgery. For Life.™