

AESTHETIC BODY CONTOURING PORTFOLIO

MICROAIRE®
For Surgery. For Life.™

PAL® Liposuction System

PAL® is a powerful, precision body contouring system designed to reduce OR time, lessen surgeon fatigue, and speed patient recovery compared to traditional liposuction.¹

System includes:

- Lightweight, autoclavable handpiece (PAL-650)
- Autoclavable cable (5006-PAL)
- Electric console (5020)

PAL-650	PAL Powered Handpiece
5006-PAL	PAL Cable
5020	PAL Electric Console



PAL® Cannulas



Extraction
for large volume
debulking procedures



Harvesting
for producing small, consistent
lobules for fat grafting



Specialty
customized for procedures
that go beyond liposuction

*Availability of disposable & multi-use cannulas within each category may vary
**Please reference the cannula brochure for all product SKU numbers

LipoFilter® Fat Processing System

- Large Volume Canister holds up to 2,500cc of lipoaspirate
- Uses gravity decantation to gently separate adipose from supernate
- Works with any liposuction system
- Sterile-packed, single-use

ASP-CAN-2S	Sterile, single-use harvesting canister. Includes filtering canister, tubing, and 60cc syringe.
ASP-CAN-NT	Sterile, single-use harvesting canister. Includes filtering canister, tubing. 60cc syringe not included.





LipoTower®

- Surgical cart for tumescent liposuction
- Adjustable flow rates of aspiration pumps with dual pumps for aspiration and infiltration
- On/Off function of suction controlled by the PAL® handpiece or a foot pedal
- Single handed adjustable I.V. pole
- Arrives fully assembled

ASP-1021 LipoTower



Contact your local MicroAire sales representative
or MicroAire customer service for more information.

1. Katz BE, Bruck MC, Coleman WP. The Benefits of Powered Liposuction Versus Traditional Liposuction: A Paired Comparison Analysis. 2001. Dermatol Surg; 27:863-967.

©2022 MicroAire Surgical Instruments, LLC
3590 Grand Forks Boulevard, Charlottesville, VA 22911
USA Telephone: 800-722-0822 FAX: 800-438-6309 EMAIL: inquiry@microaire.com
OUTSIDE USA: +434-975-8000 FAX +434-975-4134 EMAIL: inquiry@microaire.com
Follow us on Instagram and Facebook for more updates.



LIT-PAL-PORTFOLIO Rev B

MICROAIRE®
For Surgery. For Life.™